



Oxygen Powered Suction Unit with Hose

£518.74

In stock
Only **1000** left
SKU 99914

 [Add to Wish List](#)

 [Add to Compare](#)

Details

The Oxygen Powered Suction Unit with hose.

The Oxygen Powered Suction Unit is capable of fast effective suctioning of excess quantities of fluid and particles from a patient's airway. The unique on/off button allows for fast, effective suction in emergency situations. Each suction unit comes complete with a catheter ready for immediate use.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all