



Yellow Medical Waste Bin 35L with Stainless Steel Pedal

£224.40

In stock
Only **1000** left
SKU 250003

 Add to Wish List

 Add to Compare

Details

Yellow Medical Round Waste Bin (35L) with foot pedal which features a metal clamp to lock the bin lid closed.

Use:

- Clinical Waste; Pathological Waste; Infectious Waste: Waste which requires treatment or disposal by incineration in a suitably permitted or licensed facility.

Features:

- Prevents cross contamination of viral and bacterial pathogens
- Eliminates hygiene issues associated with touching the lid of a bin in areas where hygiene concerns are paramount.



Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all