



## Zoll AED Pro Defibrillator (Semi-Auto/Manual or Semi-Auto)

£3,195.73

In stock  
Only **499** left  
**SKU 400552**

Select model

-- Please Select -- ▼

 Add to Wish List

 Add to Compare

Includes:

- Soft Carry Case
- Real CPR Help to aid rescuers in delivering optimal compressions and improving the likelihood of survival
- See-Thru CPR to reduce the length of interruptions by removing compression artifact
- A three-lead cable option (cable must be ordered separately) for basic monitoring using ECG electrodes
- Unique extended USB storage that allows immediate data transfer without leaving the scene
- IP55 and a 1.5 metre drop test
- 1 Operator's Guide
- Five (5) year factory warranty

Note: Batteries and electrodes are NOT included. They must be purchased separately.



### Details

ZOLL AED Pro is designed to perform in any environment either in hospital or out in the field. The AED Pro provides the ruggedness, portability, and advanced functionality that professional rescuers and services require from an AED.

Unique to ZOLL AEDs is a real-time feedback mechanism, known as Real CPR Help, which allows rescuers to see and hear how well they are performing the rate and depth of chest compressions during CPR.

The AED Pro offers:

- Rugged construction and exceptional durability with resistance to dust and water jets (IP 55 rating)



- The only AED in the industry to withstand a 1.5-meter drop test
- Disposable or Rechargeable battery options
- Available as semi-automatic with manual override for ALS professionals
- Built-in 7-hour memory or USB memory technology
- ECG monitoring with 3-lead cable and high-quality display
- Real CPR Help with instant feedback on the rate and depth of chest compressions during CPR through ZOLL's unique one-piece CPR D-padz Electrode
- See-Thru CPR on the AED Pro allows advanced life support providers the ability to see the patient's underlying intrinsic rhythm without stopping CPR

Other important features:

- Intelligent pediatric capability with pedi-padz II Pediatric Electrodes
- Automatic self-testing to assure a constant state of readiness
- Unrivalled battery compatibility and options (long-life disposable or rechargeable) that meet your specific needs

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Frequently bought together

### You may also need





[Zoll Adult CPR-D Pads for AED Plus & Pro - Depth & Rate Sensors](#)

£132.80

[Add to cart](#)



[Zoll Soft Carry Case for AED Plus & AED Pro](#)

£141.40

[Add to cart](#)