



SURVIVAL Handy First Aid KIT

£52.15

In stock
SKU HANDY

♥ Add to Wish List
▮ Add to Compare



Details

The SURVIVAL Handy First Aid KIT is a small but serious [first aid KIT](#) -- easy to carry, pack and store – allowing you to respond to any kind of emergency.

KIT Dimensions

Height: 19cm | **Width (Closed):** 8cm | **Length (Closed):** 15cm | **Length (Open):** 34cm | **Weight:** 535g



Leading first aid design and layout

Designed to be portable enough to be with you on the go, this compact first aid kit weighs just under 500 grams, yet still contains a [SMART Snake Bandage](#) and [Pocket CPR KIT](#), as well as gloves, saline, antiseptic wipes and a resuscitation mask, shears, a first aid instruction booklet, tweezers, wound closures and lots more.

***The new and improved Handy KIT also includes a lightweight, compact emergency blanket that can be used to treat shock and hypothermia. ***

The colour-coded labelling system also enables you to clearly identify first aid contents and explains the typical use for each item. This makes it quicker and easier to respond to any emergency situation.

Easy refill and restock

The numbered labelling system tells you how many of each item you should have – alerting you to replacement requirements and making it clear when you need to refill and [restock](#).

Quality first aid product

The bag itself is coated 600D Oxford polyester, which is durable and water resistant. It's made with strong zippers and has double stitched handles. The outer KIT can be wiped down with disinfectant liquid.

***The new and improved design for the SURVIVAL Handy First Aid KIT also features the sophisticated MOLLE system for tailored attachment, new shoulder straps and redesigned zips that are sturdier than ever before, with an extra layer of material on the*



*outside that protects both zips and the KIT's contents from the elements.***

In addition to the extensive selection of first aid items, the KIT contains CPR instructions to support your first aid training in the case of any on-the-go emergency.

On-the-go first aid now more customised than ever before

Inside your KIT, you will now find a new personal space section. It's a meshed area with a zipper that allows you to include specific items or additional first aid contents suitable for your own personal needs.

Add your own insect repellent and sunscreen, or include some personal medications, blister packs or a couple of extra bandages – it's entirely up to you.

***The improved KIT is also now bigger than ever, with an extra full mesh panel, as well as two new, deep pockets. ***

Be prepared for any emergency

Ideal for adventurers, day trippers or weekend warriors, this Handy First aid KIT fits perfectly into handbags, sports bags, day packs, backpacks, rucksacks, child strollers or baby prams.

Having been on countless adventures already, it's no wonder the Handy First Aid KIT is one of Australia's leading portable first aid kits to help prepare you for any emergency.

Components - Handy First Aid KIT

- 20 x Adhesive dressings, sterile, individually wrapped
- Bandage shears, stainless steel
- Cotton gauze swabs, (3 in pack), sterile, 7.5cm x 7.5cm
- CPR Card
- First Aid Booklet
- Hypo allergenic adhesive tape, 25mm wide roll
- 2 x Nitrile gloves (1 pair in CPR KIT)
- Note pad (located at the back of the First Aid Booklet) & pencil
- 1 x SMART Bandage, heavy duty, 10cm x 2m
- 1 x Emergency Blanket
- Resuscitation face mask, disposable, incl 1 pair gloves & 1 skin cleaning wipe (in CPR KIT)
- 1 x Saline, sterile, 30ml
- 2 x Safety Pins
- 4 x Skin cleaning wipes, with cetrimide (1 with CPR KIT)
- Splinter probes, (5 per pack), sterile
- Triangular bandage, 96cm x 96cm x 136cm
- Tweezers, rust resistant with tip guard, stainless steel
- Wound closures (3 per pack), sterile

Related Products

Check items to add to the cart or select all