



CLUNGENE COVID-19 RAPID ANTIGEN NASAL SWAB SELF TEST 5PK



£8.61

Regular Price £23.98

In stock
Only **222** left
SKU CLUNGENE-
COVID-19-RAPID-
ANTIGEN

Buy 20 for £7.99 each and **save 7%**

 Add to Wish List

 Add to Compare

Included on the Australian Register of Therapeutic Goods ([ARTG 333341](#))

Very High Sensitivity

For Self-testing

Get Results at 15-20 Minutes

Pack of 5 (1 Pack = 5 Tests)



Details

CLUNGENE COVID-19 RAPID ANTIGEN NASAL SWAB SELF TEST 5/PK

The Clungene COVID-19 Antigen Rapid Test is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab or nasal swab from individuals who are suspected of COVID-19 within the first seven days of symptoms.

Included on the Australian Register of Therapeutic Goods ([ARTG 333341](#))

Very High Sensitivity

For Self-testing

Get Results at 15-20 Minutes

Pack of 5 (1 Pack = 5 Tests)



Contents

- Test Cassette
- Extraction Reagent Tube
- Sterile Swab
- Waste Bag
- Instructions for Use

Refer to Australian Register of Therapeutic Goods (ARTG [333341](#)) for full instructions. Negative test results do not exclude infection with COVID-19 (so face masks, social distancing and good hygiene practice must be maintained)



More Information

SPECIAL PRICE	£8.61
----------------------	-------

Related Products