





AEROFORM Finger Bandage Applicator Size 01

£9.11

In stock
SKU AB0120A



 Add to Wish List

 Add to Compare
Plastic Applicator



Details

Size: Size 01 **Pricing Unit:** unit/1 **Carton QTY:** 6 AEROFORM Tubular Bandage Finger Applicator is a plastic applicator used to hold the AB0120 finger tubular bandage in place, which are suited to treatment of sports injuries, dislocations and sprains.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes