



## AEROWOUND Combine Dressing 10 x 10cm Bag/50

£36.00

In stock  
SKU ACD10105

 Add to Wish List

 Add to Compare

Combine Dressing bag/50

### Details

**Size:** 10 x 10cm **Pricing Unit:** bag/50 **Carton QTY:** 16 The AeroWound Combine Dressing is recommended where a highly absorbent dressing is required for large and heavily bleeding wounds.

#### Low Adherent

- The AeroWound Combine Dressing has a special low adherent coating which minimizes trauma associated with changing wound dressings

#### Hypo-Allergenic

- The AeroWound Combine Dressing is hypo-allergenic with a low risk of wound irritation

#### Highly Absorbent

- The AeroWound Combine Dressing is capable of absorbing fluids from large and heavily bleeding wounds to stop bleeding and reduce infection risk.

#### Individually packed and sterile

- Dressing packs are available individually on demand in an easy peel sterile pouch to ensure wound cleanliness.



### Range of Sizes Available

- Ranges of dressing sizes are available to suit both small and large sized wounds effectively.





### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes