



AEROWOUND Combine Dressing 20 x 91cm

£8.63

In stock
SKU ACD20915

 Add to Wish List

 Add to Compare

Combine Dressing unit/1

Details

Size: 20 x 91cm **Pricing Unit:** unit/1 **Carton QTY:** 100 The AeroWound Combine Dressing is recommended where a highly absorbent dressing is required for large and heavily bleeding wounds.

Low Adherent

- The AeroWound Combine Dressing has a special low adherent coating which minimizes trauma associated with changing wound dressings

Hypo-Allergenic

- The AeroWound Combine Dressing is hypo-allergenic with a low risk of wound irritation

Highly Absorbent

- The AeroWound Combine Dressing is capable of absorbing fluids from large and heavily bleeding wounds to stop bleeding and reduce infection risk.

Individually packed and sterile

- Dressing packs are available individually on demand in an easy peel sterile pouch to ensure wound cleanliness.



Range of Sizes Available

- Ranges of dressing sizes are available to suit both small and large sized wounds effectively.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes