



AEROCREPE Elastic Crepe Bandage 2.5cm x 4M Wrap/12

£11.66

In stock
SKU ACE25

 Add to Wish List

 Add to Compare

Elastic Crepe Bandage 2.5cm


Details

Size: 2.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 60 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

Elastic Crepe

AEROCREPE elastic cotton crepe bandages are of an elastic conforming nature which is durable and suitable for immobilisation and compression techniques.

Easy to Open Wrapper

AEROCREPE elastic cotton crepe bandages have a easy open tab built into the wrapper of the bandage 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes