



# AERODIAGNOSTIC Personal Infrared Ear and Forehead Thermometer

£115.21

In stock  
SKU ADT05

 Add to Wish List

 Add to Compare

Ear and Forehead Infrared Thermometer



## Details

**Pricing Unit:** unit/1 **Carton QTY:** 40 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

## Infrared Thermometers

AERODiagnostic Clinical Forehead Infrared Thermometers ensure accurate and stable body temperature readings, without the need for physical contact from a distance of around 5cm.

Simply point the thermometer towards the forehead to read the temperature shown on the LCD screen. The AERODiagnostic Non-Contact Forehead IR Thermometer is specially designed to take the body temperature of a person regardless of room temperature.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes