



SUPERKID 2 Series Plastic Waterproof First Aid Kit 13 x 21 x 7.5cm

£52.95

In stock
SKU AFAK2WSK


 Add to Wish List

 Add to Compare

SUPERKID 2 Series Plastic Waterproof First Aid Kit 13 x 21 x 7.5cm



Details

Not just a kid, a Superkid! For when the invincibility powers wear off, the Superkid First Aid Kit is here! All the essentials for busy little superkids on the go! This compact kit has all the first aid items needed for treating minor injuries. Neatly packed into a tough, weatherproof case that can be wall-mounted for easy accessibility. With protection against water and dust, its a great option for use in a range of environments. Durable and well-sized for easy storage in small areas. Great for the home or car for those unexpected emergencies! 

KIT CONTENTS:

| CODE | DESCRIPTION | QUANTITY |
|-----------|--|----------|
| AGS753S-I | AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3s | 3 pc |
| AW8000-I | AEROWIPE Cleansing Wipe | 3 pc |
| AP1011-I | AEROPLAST Plastic Plaster 72x19mm | 10 pc |
| AFS001-I | AEROSHIELD Face Shield Disposable | 1 pc |
| AP325-I | AEROPORE Microporous Paper Tape 2.5cm x 5M | 1 Roll |
| APD101S-I | AEROPAD Low Adherent Dressing 10cm x 10cm | 2 pc |
| AWD14S-I | AEROWOUND Wound Dressing #14 | 1 Roll |
| AF50 | AEROFORM Conforming Bandage 5cm x 4M | 1 Roll |
| AF75 | AEROFORM Conforming Bandage 7.5cm x 4M | 1 Roll |



| | | |
|-----------|---|--------|
| ABC10 | AEROBAND Triangular Bandage 110 x 110 x 155cm | 1 pc |
| ASP01-I | AEROPINS Safety Pin | 3 pc |
| ASC13-I | AEROINSTRUMENTS Scissor 10cm | 1 pc |
| AF13-I | AEROINSTRUMENTS Tweezer 12cm Plastic | 1 pc |
| ASN096P-I | AEROSUPPLIES Notebook & Pen | 1 pc |
| AW1000-I | AEROWASH Eye Wash Ampoule 15ml | 1 pc |
| AB310-I | AEROBURN Burn Gel Sachet 3g | 1 pc |
| AGNPF02 | AEROGLOVE Nitrile Examination Gloves | 1 Pair |
| AGL001 | AEROGUIDE First Aid Leaflet | 1 pc |



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |