



DEFENDER 3 Series Metal Tough First Aid Kit 23 x 16.5 x 7cm

£114.94

In stock
SKU AFAK3M

 Add to Wish List

 Add to Compare


DEFENDER 3 Series Metal Tough First Aid Kit 23 x 16.5 x 7cm

Details

Case: Metal

Size: 23 x 17 x 70cm

Pricing Unit: unit/1

A handy kit for first response and a vital addition to any outdoor activities. Keep the whole family protected with the quality range of first aid products to treat common injuries. Ideal for travelling, the Versatile, Neat and Weatherproof cases have internal compartments to keep contents organised. For rugged conditions, opt for the Tough metal case which can also be wall-mounted. Both Tough metal and Weatherproof cases feature a rubber seal to maximise contents protection. This kit is available in 4 case options: a Tough Metal case, a Versatile Softpack reinforced fabric bag, a Weather resistant plastic case and a Neat plastic case - all our cases and bags are available with custom labelling. 

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

KIT CONTENTS:

| CODE | DESCRIPTION | QUANTITY |
|-----------|--|----------|
| AP10125-I | AEROPLAST Plastic Plaster 72x19mm | 25 pc |
| AB3-I | AEROBURN Burn Gel Sachet 3.5g | 1 pc |
| AGS753S-I | AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's | 3 pc |



| | | |
|-----------|---|---------|
| AW8000-I | AEROWIPE Cleansing Wipe | 3 pc |
| APD101S-I | AEROPAD Low Adherent Dressing 10cm x 10cm | 2 pc |
| AEP1S-I | AEROPAD Eye Pads | 1 pc |
| AF50-I | AERIFORM Conforming Bandage 5cm x 4M | 1 Roll |
| AF75-I | AERIFORM Conforming Bandage 7.5cm x 4M | 1 Roll |
| AFHP100-I | AERIFORM Heavyweight Conforming Premium 10cm x 4m | 1 Roll |
| AWD13S-I | AEROWOUND BPC Wound Dressing #13 | 1 pc |
| AWD14S-I | AEROWOUND BPC Wound Dressing #14 | 1 pc |
| AB3-I | AEROBURN Burn Gel Dressing 10x10cm | 1 pc |
| AW1000-I | AEROWASH Eye Wash Liquid Ampoule 15ml | 2 pc |
| AGNPF02-I | AEROGLOVE Nitrile Examinaion Gloves | 2 Pairs |
| AFS001-I | AEROSHIELD CPR Face Shield | 1 pc |
| ABC10-I | AEROBAND Triangular Bandage 110 x 110 x 155cm | 2 pc |
| AP325-I | AEROPORE Microporous Tape 2.5cm x 5M | 1 Roll |
| ASC13-I | AEROINSTRUMENTS Scissor 13cm | 1 pc |
| AF11-I | AEROINSTRUMENTS Tweezer 11cm | 1 pc |
| ATB130-I | AERORESCUE Emergency Rescue Blanket | 1 pc |
| ASN096P-I | AEROSUPPLIES Notebook & Pen | 1 pc |
| ASP12-I | AEROPINS Safety Pin | 3 pc |
| AGL001-I | AEROGUIDE First Aid Leaflet | 1 pc |



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |