



RESPONDER 4 Series Plastic Waterproof First Aid Kit 24.3 x 24.3 x 8cm

£164.91

In stock
SKU AFAK4W

 Add to Wish List

 Add to Compare


RESPONDER 4 Series Plastic Waterproof First Aid Kit 24.3 x 24.3 x 8cm

Details

Case: Plastic

Size: 26 x 36 x 9cm

Pricing Unit: unit/1

Ensure workplace compliance for first aid with this comprehensive range, containing a wide scope of products to cater for any workplace emergencies. Respond promptly to incidents with a variety of outer cases to suit all types of environments. Internal compartments and shelves keep product organised, and all plastic and metal case options are able to be wall mounted. The softpack style is ideal for mobile workers, while the durable plastic cases are recommended for industrial areas. **Contents meet Australian Standard AS2675-1983**  This kit is available in 5 case options: a Robust rugged plastic case, Versatile Softpack reinforced fabric bag, a Weather resistant plastic case and a Neat plastic case - all our cases and bags are available with custom labelling.

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

KIT CONTENTS:

CODE	DESCRIPTION	QUANTITY
AP10125-I	AEROPLAST Plastic Plaster 72x19mm	25 pc



AFP50125-I	AEROPLAST Fabric Premium Strip 72mm x 19mm	25 pc
AB3-I	AEROBURN Burn Gel Sachet 3.5g	5 pc
AAS50-I	AEROAID Antiseptic Spray 50ml	1 pc
AFS001-I	AEROSHIELD CPR Face Shield	1 pc
AGNPF02-I	AEROGLOVE Nitrile Examination Gloves	5 Pairs
AGS753S-I	AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's	5 Pks
AW1000-I	AEROWASH Eye Wash Ampoule 15ml	8 pc
AW8001-I	AEROWIPE Cleansing Wipe	10 pc
ASP001-I	AEROPROBE Splinter Probes 3.7cm	10 pc
APD50S-I	AEROPAD Low Adherent Dressing 5cm x 5cm	6 pc
APD100S-I	AEROPAD Low Adherent Dressing 7.5cm x 10cm	3 pc
APD101S-I	AEROPAD Low Adherent Dressing 10cm x 10cm	2 pc
AF50-I	AEROFORM Conforming Bandage 5cm x 4M	3 Roll
AF75-I	AEROFORM Conforming Bandage 7.5cm x 4M	3 Roll
AFHP100-I	AEROFORM Heavyweight Conforming Bandage 10cm x 4m	1 Roll
AP325-I	AEROPORE Microporous Tape 2.5cm x 5M	1 Roll
AWD14S-I	AEROWOUND BPC Wound Dressing #14	1 Roll
AWD15S-I	AEROWOUND BPC Wound Dressing #15	1 Roll
ACD1020S-I	AEROWOUND Combine Dressing 10cm x 20cm	1 pc
ABC10-I	AEROBAND Triangular Bandage 110 x 110 x 155cm	2 pc
ATB130-I	AERORESCUE Emergency Rescue Blanket	1 pc
AEP1S-I	AEROPAD Eye Pads	4 pc
AI2500-I	AEROCOOL Instant Ice Pack 80g	1 pc
ASC13-I	AEROINSTRUMENTS Scissor 13cm	1 pc
AF13-I	AEROINSTRUMENTS Tweezer 11cm	1 pc
AAB100-I	AEROSUPPLIES Amputated Parts Bags Assorted Sizes	3 pc
ASN096P-I	AEROSUPPLIES Notebook & Pen	1 pc
ASP12-I	AEROPINS Safety Pin	12 pc
AGL001-I	AEROGUIDE First Aid Leaflet	1 pc



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes