




MODULATOR Blue Eye Wound Module 20 x 10 x 6cm

£37.92

In stock
SKU AFAKMODE

 Add to Wish List

 Add to Compare

MODULATOR Blue Eye Wound Module 20 x 10 x 6cm



Details

The Modulator Eye Wound Module is a visual and clearly labelled soft pack that can be ordered individually as a replacement module for The Modulator Softpack or The Modulator Metal Cabinet editions.

This module can also be used in addition to any core workplace kit for compliance where a high risk of eye injury has been identified or as a single standing eye wound kit.

Ensures hasslefree servicing and easy re-stocking with individual replacement modules.

Contents meet Australian Standard AS2675-1983



CODE	DESCRIPTION	QUANTITY
AP312	AEROPORE Microporous Paper Tape 1.25cm x 5M	1 Roll
AF75	AEROFORM Conforming Bandage 7.5cm x 4M	1 Roll
AEP1S	AEROPAD Eye Pads 5.5cm x 7.7cm	6 pc



AW1000	AEROWASH Eye Wash Ampoule 15ml	10 pc
AGEL	AEROGUIDE Eye Wound Treatment Card	1 pc
AGNPF02	AEROGLOVE Nitrile Examination Gloves Pair	1 pc



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes