



AEROFORM Heavyweight Conforming Bandage 7.5cm x 4M Wrap/12

£33.17

In stock
SKU AFH75

 Add to Wish List

 Add to Compare

Heavy Weight Conforming Bandage 7.5cm

Details

Size: 7.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 30 The AeroForm Conforming bandage is a quality cost-effective bandage with excellent conforming attributes. Our Heavyweight and Premium range offers greater support and compression and are recommended for sprains, strains and snake bites.

Conforming Specifications

AeroForm Heavy Weight Conforming Bandages are of a heavy weight conforming nature which is used to support sprains and strains and to provide comfort to injured limbs and joints.

Easy to Open Wrapper

AeroForm Conforming Bandages have a easy open tab built into the wrapper of the bandage. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes