



## AEROFIX Non-Woven Island Dressing 9 x 15cm Box/3

£8.21

In stock  
SKU AFID153S

 Add to Wish List

 Add to Compare

Aerofix non-woven Island Dressing



### Details

**Size:** 9 x 15cm **Pricing Unit:** env/3 **Carton QTY:** 96 AEROFIX dressings combine a SupaSorb pad with a soft and conformable fixative layer for the simple and effective management of larger wounds, abrasions and sutured wounds.

#### SupaSorb Wound Pad

- An embossed wound pad provides superior wound management by draining blood or exudates into pockets in the wound pad. The Super absorbent wound pad minimises the number of dressing changes whilst protecting the wound from further injury. SupaSorb™ leaves the wound site clean and comfortable, while the non-adherent film reduces the pain and site trauma when changing the dressing.

#### SupaBond adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

#### Soft and Conformable

- A non-woven backing is smooth and soft, and adheres to body contours using SupaBond adhesion.



### Patient Comfort

- The SupaSorb non-adherent wound pad , combined with a soft and conforming non-woven backing takes patient comfort to a new level, minimising trauma and assisting recovery.

### Contraindications / Precautions

- AEROFIX should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This product is not a substitute for direct medical attention from a healthcare practitioner.** 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes