



## AEROFIX Underwrap Tape 5cm x 10M

£6.76

In stock  
Only **671** left  
**SKU** AFR50



 Add to Wish List

 Add to Compare

Aerofix Non-Woven Fixation Tape



### Details

**Size:** 5cm x 10M **Pricing Unit:** unit/1 **Carton QTY:** 100 AeroFix non-woven fixation tape is a highly conformable, latex-free solution for dressing retention and underwrap.

#### SupaBond adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

#### Soft and Conformable

- A non-woven backing is smooth and soft, and adheres to body contours using SupaBond adhesion.

#### Contraindications / Precautions

- AeroFix should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This product is not a substitute for direct medical attention from a healthcare practitioner.** 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes