



## AEROTAPE White Microporous Paper Tape 2.5cm x 5M Box/12

£7.97

In stock  
SKU AP325



 Add to Wish List

 Add to Compare

Microporous Paper Tape



### Details

**Size:** 2.5cm x 5M **Pricing Unit:** box/12 **Carton QTY:** 30 AEROPORE Microporous Paper Tape is a flexible and breathable paper tape suitable for securing dressings, especially to facial and fragile skin. AeroPore Non-woven Microporous Paper Tape is recommended for the securing of dressings and medical equipment on sensitive and fragile skin including facial areas.

### Hypoallergenic

- AEROPORE Microporous Paper Tape has a very low irritant risk, and is suitable for use on fragile and delicate skin.

### Individually Wrapped

- AEROPORE Microporous Paper Tape is available in boxes of individually wrapped rolls ensuring the tape is kept in optimum condition.





### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes