



## Green Terylene/Cotton Carry Sheet



£53.14

In stock  
**SKU** AR1017

 Add to Wish List

 Add to Compare

Green Cotton Carry Sheet Terylene



### Details

**Size:** 190 x 68cm **Pricing Unit:** unit/1 **Carton QTY:** 1 The green cotton terylene carry sheet is a compact and convenient carry sheet suitable for a range of actions for fast patient recovery. This product is highly suitable for public emergency with its key attributes being flexibility and strength. The carry sheet folds to a small package; yet has 8 strong carry handles for patient handling.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes