



AEROSPLINT Rolled Aluminium Foam Splint 90 x 11cm

£14.05

In stock
SKU AR1018R

 [Add to Wish List](#)

 [Add to Compare](#)

Aluminium Foam Splint Rolled

Details

Size: 90 x 11cm **Pricing Unit:** unit/1 **Carton QTY:** 50

Light Weight Rolled Splint

Immobilisation is a critical item of care when responding to a casualty to limit the occurrence of further injury. Quality immobilisation using suitable products is integral to a full recovery of skeletal injuries.

Aluminium Foam Splints

An aluminium foam splint is a highly adaptive splint that can be rolled folded or formed into almost any application with the use of a quality support bandage. Available in a rolled or flat form, shrink-wrapped for protection.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes