




AEROWOUND Finger Dressing Bag/50

£35.70

In stock
SKU AWF002

 Add to Wish List

 Add to Compare
Finger Dressing



Details

Size: 4.5 x 4.5cm **Pricing Unit:** bag/50 **Carton QTY:** 20 The AEROWOUND Wound Dressings are a general purpose wound dressing for use in preventing infection and controlling bleeding. Available in sizes to suit wounds from finger injuries to large wounds.

Highly Absorbent

- Highly absorbent wound pad controls bleeding and absorbs wound exudates

Low adherent

- Special low adherent coating minimizes wound trauma associated with changing dressings

Individually packed and sterile

- Dressings are available individually on demand in an easy tear wrapper and sterile for wound cleanliness.

Ease of application

- The wound pad is attached to a wrapping leader of gauze bandage making the dressing simple to apply, yet comfortable and secure for all general wound dressing.



Unique Finger Dressing

- The AEROWOUND Finger dressing is unique in the industry for its ability to treat finger wounds too large to cover with an adhesive dressing. This unique product combines a 4.5cm square low adherent and absorbent pad, with a wrapping leader and secure adhesive tab to protect the wound.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes