



CARDIACT Connect Monitored AED Cabinet

£1,630.45

In stock
SKU CC-300

 Add to Wish List

 Add to Compare

Outdoor Wall Cabinet with Monitoring for any AED



Details

Size: -

Pricing Unit: unit/1

Carton QTY: 1



Outdoor Wall Cabinet with Monitoring for any AED



[dfliip id="11144" type="thumb"][/dfliip]



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes