




## ChoiceMMed Nylon Carry Case for Pulse Oximeter



£6.55

In stock  
**SKU** CNC

 Add to Wish List

 Add to Compare

Nylon Carry Case for Pulse Oximeter



### Details

**Size:** Fits C63 and C22 **Pricing Unit:** unit/1 **Carton QTY:** 1 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

### Soft Design

Designed and built with a soft nylon finish.

### Durable

Built soft but with a very thick and durable nylon and will keep your oximeter safe

### Compatability

Compatible and designed for most oximeters

### Quality

Built with great stitching and a strong nylon material.



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes