



LIFEPAK CR2 Fully-Automatic Defibrillator with Wi-Fi

£2,210.38

In stock
SKU CR-2-A

 Add to Wish List

 Add to Compare

Fully-Automatic Defibrillator AED WiFi

Details

Size: - **Pricing Unit:** unit/1 **Carton QTY:** 1

- 8-year warranty
- Wireless connectivity
- Fastest time to 1st shock
- CPR Insight
- Durable
- 4 years of standby battery life

Also available as a [Semi-Automatic](#) AED. Sudden cardiac arrest (SCA) can happen to anyone, anywhere. Immediate treatment is vital. A victim's chance of survival dramatically decreases for every minute without treatment. That's why public access defibrillators are so important: they put lifesaving technology where it can do the most good. So when an emergency happens, you should have nothing less than the best. In a cardiac emergency, the LIFEPAK CR2 defibrillator connects everything and everyone involved in the response reducing delays in treatment and improving care across the lifesaving system.


Who will be saved next?

When you choose the LIFEPAK CR2 defibrillator, you're choosing a connected solution that empowers bystanders to help save a life all while backed by innovative technology that allows continuous CPR and streamlines communication across



the system of care.

With a design influenced by user feedback, it's no surprise the CR2 is rated highest in overall user confidence by AED users. The groundbreaking LIFEPAK CR2 is at the heart of a complete AED response system.

 **Layered design** Layered design with easy-to-follow bold graphics. Both trained and untrained AED users clearly know how to begin. **QUIK-STEP™ electrodes** Peel directly off the base for faster placement. **cprINSIGHT analysis technology** Analyses for shockable rhythm during chest compressions with no need to pause. **Metronome and CPR coaching** Sets an effective pace and audibly guides users. **ClearVoice™ technology** Detects background noise and adjusts volume accordingly for clear instruction. **Fully automatic** Available in fully or semi-automatic models. **Highest available escalating energy** Up to 360J for more effective shocks as needed. **Bilingual** Toggle between two pre-set languages when using the device. **Child mode** Child mode delivers lower energy levels appropriate for young children without having to change electrodes. **LIFEPAK TOUGH™** IP55 rating for challenging environments. **8-year warranty** Backed by an 8-year warranty. **LIFELINKcentral AED Program Manager** Monitor AED programs by tracking AED status, sending patient data to emergency responders and hospitals, detecting AED locations and other tools. [Register your Product Warranty Here.](#)



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes