



Disposable Green Plastic Kidney Dish 200mL




£1.61

In stock
SKU KD150D

Buy 100 for £1.35 each and **save 16%**
Buy 200 for £1.30 each and **save 19%**
Buy 1000 for £1.11 each and **save 31%**

 [Add to Wish List](#)

 [Add to Compare](#)
Plastic Kidney Dish



Details

Size: 200mL **Pricing Unit:** unit/1 **Carton QTY:** 100 The AEROSUPPLIES stationary range covers a number of common stationary requirements applicable to the First Aid Industry.

Disposable Plastic Kidney Dish

Kidney dishes are a widely used product in surgical procedures for holding small instruments and more importantly for rinsing with sterile water. The dishes are also used to contain samples for laboratory analysis. The AeroSupplies Plastic Kidney dish is disposable. We are able to source autoclavable kidney dishes if you have a specific request for reusable kidney dishes.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes