



AERORESCUE Alloy Dual-Fold Emergency Pole Stretcher

£173.34

In stock
SKU STR-02



 Add to Wish List

 Add to Compare

Alloy Emergency Pole Stretcher



Details

Size: 207 x 54 x 14cm **Pricing Unit:** unit/1 **Carton QTY:** 1 AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-02

Weight limit: 160kg Popular basic stretcher offering great value for money. Folds to 103cm x 18cm x 10cm in a convenient carry bag



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes