







AERORESCUE Alloy Scoop Stretcher

£462.75

In stock
SKU STR-06



 Add to Wish List
 Add to Compare
Alloy Scoop Stretcher



Details

Size: 190 x 43 x 6cm **Pricing Unit:** unit/1 **Carton QTY:** 1 AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-06

Weight limit: 160kg Scoop stretcher designed for gently scooping up a patient while maintaining the found position minimising the risk of further injury. Extendable head rest for 170cm-190cm patient size. Folds down to 172cm x 43cm x 6cm



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes