



QV Intensive Cream - 500g

£14.59

In stock
Only **498** left
SKU 11372

 [Add to Wish List](#)

 [Add to Compare](#)

QV Intensive Cream is a rich cream scientifically designed to moisturise and help maintain moisture in extremely dry skin. The preservative-free formulation is suitable for very sensitive skin, or skin prone to sensitivity due to preservatives.



Details

QV Intensive Cream is a rich cream scientifically designed to moisturise and help maintain moisture in extremely dry skin. The preservative-free formulation is suitable for very sensitive skin, or skin prone to sensitivity due to preservatives.

Benefits:

- Triple-moisturising action, including 10% glycerin, for optimal skin hydration without a greasy residue
- Clinically proven 24-hour moisturisation
- Easily applied and readily absorbed into the skin without excessive rubbing
- Helps maintain the skin's natural protective barrier and lock in moisture
- Preservative-free, further reducing the chance of irritancy



- Low-irritant formulation, fragrance, lanolin, propylene glycol and colour-free
- Non-comedogenic, so it won't block pores
- Hypoallergenic & Dermatologically tested
- Vegan Friendly - contains no animal-derived ingredients
- Cream suitable for use with sensitive and eczema-prone skin

Directions:

Apply Intensive Cream to affected areas as often as required, particularly after washing, showering, bathing, shaving, exposure to harsh climatic conditions and at night.

Cautions:

Avoid eyes. For external use only. Store below 30°C.

Ingredients:

Aqua (Water), Petrolatum, Paraffinum Liquidum, Glycerin, Cyclopentasiloxane, Cyclohexasiloxane, Cetearyl Alcohol, Behentrimonium Methosulfate, Dimethicone, 1,2-Hexanediol, Caprylyl Glycol.