



## Clinell Universal Sanitising Wipes with Clip Packet 50

£6.47

In stock  
Only **100** left  
**SKU** CWCP50AUS



 Add to Wish List  
 Add to Compare



### Details

## Clinell Universal Wipes with Clip PKT/50

**Choose Clinell Universal Wipes as the wipe that cleans and disinfects in one easy step.**

Clinell Universal Wipes are a combined cleaning and disinfectant wipe, suitable for use on surfaces, killing at least 99.99% of germs, including Coronavirus.

Perfect for use on surfaces and equipment that are frequently used including tabletops, door handles, light switches, lift buttons, desks, keyboards, chairs, check-outs, EFTPOS machines, bathrooms, and kitchens.

- Clinell Thick Wipes - extra detergent and thicker wipes
- Effective against bacteria and viruses, including COVID-19
- Hospital Grade Disinfectant
- Multipurpose cleaning and disinfecting solution
- Kills at least 99.99% of bacteria and viruses
- Dermatologically tested
- No bleach, no taint
- Safe for kitchen areas
- ARTG listed 178363



- **Unit** PKT/50

## **Clinell Universal Wipes are effective against COVID-19, MRSA and VRE.**

### **More information:**

- [Clinell universal TGA.PDF](#)
- [Clinell Brochure.pdf](#)
- [Cleaning\\_and\\_disinfecting\\_during\\_Covid-19.pdf](#)

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare providers. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare providers before starting any new treatment or with any questions that you may have regarding a medical condition or injury.