



## La Petite Spacer

£11.65

In stock  
Only **49** left  
**SKU** BHC-LP-1

 Add to Wish List

 Add to Compare

The e-chamber La Petite is a spacer for effective delivery of medication suitable for all ages. The compact 220ml separable chamber is made from durable polypropylene plastic, which is entirely recyclable. This is also anti-static to ensure minimum medication adheres to the chamber, unlike alternatives.



### Details

Compact asthma spacer ideal for people of all ages.

Features include:

- 220ml
- Anti-static plastic with embedded particles - can't be scratched or washed off
- Single patient re-usable
- Easy cleaning
- Universal inhaler fitting
- 15mm-22mm round mouthpiece
- 100% recyclable
- Inhaler can be conveniently stored inside
- Retail packed and available in quantities of 1, 10, 25 or 50 per box