




e-Chamber Peak Flow Meter

£15.37

In stock
Only **30** left
SKU BHC-PFM

 Add to Wish List

 Add to Compare

CLASSIC FULL RANGE PEAK FLOW METER

The e-chamber Peak flow meter is a personal healthcare device designed for home respiratory monitoring with full lung volume range (60-900L/min)



Details

Technical Specification:

- Measuring Range: 60-900 L/MIN (EU & ATS)
- Scale increment: 10L/MIN (over whole scale)
- Calibration: 100% (individual instrument)
- Accuracy: $\pm 5\%$ or 10L/MIN Inter-instrument
- Repeatability: $\pm 2\%$ Intra-instrument
- Flow Impedance: 0.30kPa/l/s@c550l/min
- Dynamic State: Critically damped @ 350L/MIN
- Frequency Response: 100 Hz @ 300L/MIN
- Standard Compliance: ISO 23747; EN 13826; AS/NZ 4237; NAEP Guide 1991 (ATS)
- EU Directive: 93/42/EEC