



Medium Sports Tape

£8.39

In stock
Only **998** left
SKU 101458

 Add to Wish List

 Add to Compare

Medium Sports Tape 38mm x 15mm

Product Details:

- Brand: Trafalgar
- Colour: Beige
- Size: Medium

Features:

- Hand Tearable
- Zinc Oxide Ridged Tape
- Strong Support
- Suitable For Shoulder, Wrist, Knee, Ankle And Foot

Details

Trafalgar Medium Sports Tape

Comfort and safety when you are out and about

Trafalgar's sports tape is made from strong woven fabric that has zinc oxide adhesive for strong adhesion.

The Sports Tape comes with serrated edges for easy tearing - meaning that you do not have to have scissors on hand!

This strapping tape is highly durable and helps to support and restrict movement in joint areas.

The breathability of this tape means that sweat can evaporate, maintaining the strength and stability of the tape whilst you are doing physical activity or sport.

Medium Sports Tape 38mm x 15m