



## Cohesive Bandages Tan 5cm

£6.47

In stock  
Only **996** left  
**SKU** cohesive-  
tan-5cm



 Add to Wish List

 Add to Compare

- Great for the kitchen
- Use as a strapping tape for sport
- Reusable



### Details

## Cohesive Bandage Tan 5cm

*Maintains compression and is great for sports strapping*

- Great for the kitchen
- Use as a strapping tape for sport
- Tan coloured to be more discrete
- Reusable

Cohesive bandages are unique in that they do not need any adhesive, clips or pins to secure. This means ease when applying and taking off, as these bandages will not stick to skin or hair.

These hand-tear Tan coloured bandages can be carried with you, or be a useful addition to your first aid kit without the need for any additional cutting materials like scissors.



With special lightweight and breathable technology, this bandage can act as a secondary dressing.

Cohesive bandages are unique in that they maintain the compression/tension that you choose.

Cohesive bandages with their tension retention can offer support for weak or injured joints e.g. sprained ankle, knee ligament injury

Unlike other strapping tapes, these cohesive bandages maintain their tension meaning that you do not have to make any adjustments when you apply

- Applying the bandage with stretch gives extra tension
- Applying the bandage without stretch provides light support
- Extra layers can of bandage can provide more support

Dimensions: 5cm x 4.5m

**Contents:**

1 x roll



**More Information**

**SPECIAL PRICE**

£5.66