



## Cotton Applicator 3's (100) 7.5cm Sterile Single Tip

£50.02

In stock  
Only **999** left  
**SKU** 7515

 Add to Wish List  
 Add to Compare

### Details

Cotton Tips are soft, absorbent and hygienic for medical and general purpose use.

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.