





## Melolite Low Adherent Dressing (100) 7.5 x 20cm - Low

**\$95.27**

In stock  
Only **1000** left  
**SKU** 4942

 Add to Wish List  
 Add to Compare

### Details

Melolite dressings are constructed from an absorbent pad which helps distribute the fluid, thus avoiding the problems of immediate exudate strike-through, covered on both sides with a heat-laminated polyester film. This film acts as a low-adherent wound contact layer by allowing the rapid passage of exudate through its elliptical holed construction into the absorbent pad. This minimises the risk of adherence to the wound and the related trauma on dressing removal. Either side of the dressing can be placed directly onto the wound site.

### Features and Benefits:

- Double sided low-adherent wound dressing
- Absorbent pad
- Highly conformable to all wound sites
- Sterile Presentation

**Quantity:** 100

**Dimensions:** 7.5cm x 20cm

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



---

**Related Products**

Check items to add to the cart or select all