





## Najo Redihold Backboard with Pins - 1830 x 410 x 40mm

**\$621.74**

In stock  
Only **1000** left  
**SKU** 280118

 Add to Wish List  
 Add to Compare

### Details

The Najo RediHold is our top of the line backboard with two handles at a tapered end for stable transport. The head area is smooth to insure all head immobilisation devices adhere to the board and the edges are angled for log-rolling.

Strong, lightweight, one-piece HDPE construction is impervious to all bodily fluids, is easy to decontaminate, and provides 100% X-ray translucency. This board is also buoyant and can be used to float a patient during water rescue operations.

Guaranteed for life under normal usage.

- Dimensions: 183 x 41 x 4cm
- Weight: 4kgs
- Load Limited: 272kgs

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all