



Rediwipe Isopropyl Wipes Cannister (75) 14.5cm x 42cm

\$20.12

In stock
Only **140** left
SKU 99449

 Add to Wish List
 Add to Compare



Details

Rediwipe Isopropyl Wipes Cannister (75) 14.5cm x 42cm

- 70% Isopropyl alcohol impregnated wipes
- 75 wipe per cannister
- Multi purpose wipe
- Designed for sanitising clean surfaces
- Excels in removing adhesive and gums

Dimensions: 14.5cm x 42cm

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products



Check items to add to the cart or select all