





Salinaax Emergency Eye & Skin Wash Station with Saline - 2 x 944ml Saline soft bottles

\$58.40

Out of stock
Only **1000** left
SKU 260011

 [Add to Wish List](#)
 [Add to Compare](#)

Details

The Eye Wash Station is designed as a secondary eye wash device and is intended for support use, not as a replacement for self contained or plumbed units.

The Eye Wash Station should be located in the immediate vicinity of employees working in a potentially hazardous area. It is used to deliver initial first aid while getting the person to a primary unit, or as extended aid when transporting a person from the primary eye wash unit to a medical facility.

Time is critical, the range of eye wash bottles have been designed to ensure that opening and commencement of flushing is achieved in optimal time. Just twist the top to break the seal and commence flushing.

Contents:

- Salinaax Emergency eye and skin wash station
- 2 x 944ml saline soft bottles

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



Related Products

Check items to add to the cart or select all