



## Sam Finger Splint Universal 9.5 x 4cm

**\$3.18**

In stock  
Only **1000** left  
**SKU** 490003

 Add to Wish List  
 Add to Compare

### Details

The Sam Splint Original is an extremely pliable multi-purpose splint and is built from a thin core of aluminium alloy, sandwiched between two layers of closed cell foam. Bent into various shapes and curves, it becomes extremely strong and supportive for any fractured or injured limb. The Sam Splint Original is also ideal for splinting IV's.

It is lightweight, radiolucent, reusable and very portable. Almost any bone in the body can be splinted with Sam Splint.

- Size: 9.5 x 4cm

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all