



## Standard Skeleton on Stand (Model A10)

**\$476.89**

In stock  
Only **1000** left  
**SKU 8389**

 Add to Wish List

 Add to Compare



### Details

Standard Skeleton on Stand (Model A10).

The standard model of a human skeleton, has been appreciated throughout the world for decades.

This German quality model is hand assembled for durability and is one of the most widely sold skeletons across the globe. The arms and legs are removable for individual study and the skull can be dissected into three parts.

It is made of durable, washable, unbreakable plastic. Life-size human skeleton model for anatomical study. Its very high quality and robust construction.

The other advantages of the skeleton are:



- Exceptional value for money
- 3 year guarantee
- Top quality natural casting
- Final assembly carried out by hand
- Close to the realistic weight of around 200 bones

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all