



## Stax Finger Splint - Durable Plastic

**\$9.90**

In stock  
Only **1000** left  
**SKU** WEB-DSFS

 [Add to Wish List](#)

 [Add to Compare](#)



### Details

- Protects distal phalanx and extends the DIP joint.
- Open tip permits examination of fingertip and nailbed.
- Made of perforated, translucent semi-hard plastic.
- Latex free.
- Supports the DIP joint in neutral to help correct mallet deformity.
- Splint can be trimmed with heavy duty scissors.

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all