



## Stethoscope Arden Double Head

**\$42.49**

In stock  
Only **1000** left  
**SKU** 99651

 Add to Wish List  
 Add to Compare



### Details

Stethoscope Arden Double Head is a lightweight, nurses style stethoscope with bell and diaphragm. This stethoscope is good for a basic general exam listening to breathing sounds, bowel sounds, and taking blood pressure.

### Specifications:

- Tube length: 25inches
- Bell diameter: 3.3cm
- Diaphragm diameter: 4.4cm

### Options:

- Black
- Blue
- Red

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical



---

emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

**Related Products**

Check items to add to the cart or select all