





## Tegaderm Transparent Film Dressing Plus Pad 5cm x 7cm

**\$79.25**

In stock  
Only **1000** left  
**SKU 230013**

 [Add to Wish List](#)  
 [Add to Compare](#)



### Details

This composite dressing is a sterile, waterproof bacterial barrier which consists of a non-adherent absorbent pad bonded to a larger thin film transparent dressing. Tegaderm +Pad transparent dressings, with the unique "picture frame" delivery system, make it quick and easy to achieve a precise, secure placement of the dressing every time.

**Size:** 5cm x 7cm

**Quantity:** 50

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.