





## Thermometer Infra-red Ear

**\$111.12**

In stock  
**SKU** 99553

 Add to Wish List  
 Add to Compare

### Details

Bodicherk Thermometer Infra-red EarSafe is accurate and easy to use.

#### Features:

- Fast readings in 0.1 seconds
- Infrared technology for accurate readings
- Easy one touch operation
- Memory recalls last 10 temperature readings
- Automatic shut-off
- Convenient carry case

The Infrared Ear Thermometer utilises the latest technologies to measure infrared energy emitted from the eardrum. The ear is a highly accurate site for measuring the core body's temperature because the eardrum shares blood vessels with the hypothalamus which is the area of the brain that regulates body temperature.



---

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all