



Versa Towel Small 24.5 x  
41.5cm (100)

\$234.47

Out of stock  
SKU 4210

 Add to Wish List

 Add to Compare



## Details

- High quality, crepe wadding provides excellent absorbency
- Gentle on skin when used for repeated hand drying without chafing
- Easy and economical perforated roll format
- Unique product for medical areas that demand excellence
- Measurement: 24.5cm x 41.5cm
- Carton buy of 16
- Box of 100

## Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

## Related Products

Check items to add to the cart or select all