



## Wheelchair Lightweight Transit - Max Load 100kg

**\$288.31**

Out of stock  
**SKU** 250062

 Add to Wish List

 Add to Compare

### Details

This lightweight transit duty wheelchair features the following:

- Lightweight design
- Stainless steel frame
- Removable foot plates
- Removable arm rests
- Dual side brakes
- Standard mag rear wheel with pneumatic tyre
- Folding frame for convenient transport
- Synthetic, washable padded upholstery
- Tipping levers for the attendant
- 100kgs weight capacity

### Specifications:

- Seat width: 46cm
- Seat depth: 40cm



- Back Style: Fixed
- Arm Styles: Removable
- Front Riggings: Swing-away
- Rear Axle: Fixed position
- Rear Wheel type: Composite mag
- Frame type: Stainless Steel (powder coated)
- Upholstery: Vinyl
- Product Weight: 16kgs
- Product Width: 62cm
- Overall Length: 62cm
- Overall Height: 90cm
- Weight Capacity: 100kgs

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all