



S+M Tourniquet Latex Free, Single Use

\$6.62

In stock
Only **999** left
SKU
AXP10180001



 Add to Wish List

 Add to Compare

- Single patient use
- One-hand operation
- No pinching or jamming
- Instant removal with the quick release buckle
- Dispenser box for easy access
- Individually wrapped
- Latex free



Details

S+M Tourniquet is a latex-free tourniquet designed to be used when taking blood samples. S+M Tourniquet features a plastic clip that prevents any pinching or jamming when used on patients.

Quick and easy to use, S+M Tourniquet is soft and comfortable thus ensuring patient comfort.

- Single patient use
- One-hand operation
- No pinching or jamming
- Instant removal with the quick release buckle
- Dispenser box for easy access
- Individually wrapped
- Latex-free



Note: Price is per individual tourniquet.