



RapidStop Tourniquet

As low as **\$61.75**

In stock

[Only 200 left](#)

Product
availability

Product Name	Qty
RapidStop Tourniquet-Red	100.0000
RapidStop Tourniquet-Black	100.0000

SKU RS-
TOURNIQUET

configurable phtml

 Add to Wish List

 Add to Compare

The RapidStop® Tourniquet relies on key technologies to achieve rapid, intuitive, and easy one-hand self-application. Application is even faster when a responder uses both hands to apply the tourniquet to a victim.



Details

Intuitive

The RapidStop® Tourniquet relies on key technologies to achieve rapid, intuitive and easy one-hand self-application. Application is even faster when a responder uses both hands to apply the tourniquet to a victim.

Strong

The RapidStop® Tourniquet achieves full occlusion reliably and rapidly. It requires only gross motor control, allowing self-application to be readily possible. The RapidStop® Tourniquet is appropriate for controlling hemorrhage for injuries to both the arms and legs.

Lifesaving

Easy as 1-2-3 Three simple and fast steps set RapidStop® Tourniquet apart from the competition. Guided by instructions on the product, pocket-sized instruction cards and video format.

Integrated instructions on the RapidStop® Tourniquet are included using stitched waterproof vinyl tags resistant to any wear and tear.



**For use as a life-saving device in high-stress trauma situations where application time is critical.
Only gross motor skills required to operate the RapidStop® Tourniquet.**

- Police, firefighters and EMTs are often first on the scene in civilian medical emergencies. RapidStop® Tourniquet can be easily carried and deployed by these service personnel to save lives.
- In places where people work daily with dangerous equipment on a construction site, oil rig, lumber yard, machine shop or distribution warehouse, RapidStop® Tourniquet should be at hand to mitigate the risk of death through serious limb injury.
- Whether hunting, camping, rock climbing, or mountain biking for leisure, the risks of extreme trauma due to accidents are real. RapidStop® Tourniquet should be an essential piece of kit for any outdoor adventurer.



More Information

BRAND	RapidStop
-------	-----------