



## EgoPsoryl TA 30g

**\$12.89**

In stock  
Only **461** left  
**SKU 11515**

 [Add to Wish List](#)

 [Add to Compare](#)

EgoPsoryl TA is a soothing corticosteroid-free gel, formulated to help control stable, moderate psoriasis and persistent dermatitis.

### Details

EgoPsoryl TA is a soothing corticosteroid-free gel, formulated to help control stable, moderate psoriasis and persistent dermatitis.

### Benefits

- Aids in the control of stable, moderate psoriasis (including psoriasis of the scalp) and persistent dermatitis.
- Soothing gel formulation calms the skin and relieves itching.
- Non staining.

### Directions

Apply to affected skin 2 – 4 times daily. Avoid soap as soap may irritate skin. Cleanse with Pinetarsol, an anti-itch cleansing treatment which helps reduce inflammation and remove scales. For enhanced effect, apply petroleum jelly over EgoPsoryl TA.

Scalp: Apply at night if desired leave on during the day.

Cleanse with 'Sebitar' when needed.

Sensitive skin (eg face, children): Use Egoderm Cream for first 7 days. Apply EgoPsoryl TA over Egoderm Cream for 3 days. Then use EgoPsoryl TA alone.

### Ingredients

Coal Tar Solution 5%, Phenol 0.5% and Sulfur 0.5%

**Contains alcohol.**