



QV Baby Moisturising Cream Pump - 500g

\$21.64

In stock
Only **497** left
SKU 11674

 Add to Wish List

 Add to Compare

Replenishes the skin and helps relieve dry areas



Details

QV Baby Moisturising Cream is suitable for use every day and with conditions such as eczema, dermatitis and psoriasis. Contains Squalane, a component found in the skin's natural oil, helps seal in the skin's own moisture.

Directions:

Apply to the face and body, massage in gently. For best results combine with use of QV Baby Bath Oil and QV Baby Gentle Wash. Tip - Apply to affected areas as often as required, particularly after washing, showering, bathing, exposure to harsh climatic conditions and at night. For the best results combine with use of QV Baby Bath Oil and QV Baby Gentle Wash.

Caution:

Avoid eyes. For external use only. If irritation occurs, discontinue use. Store Below 30°C.



Ingredients:

Aqua (Water), Glycerin, Paraffinum Liquidum, Petrolatum, Cetearyl Alcohol, Dimethicone, Niacinamide, Stearic Acid, Cetareth-20, Squalane, Laureth-3, Glyceryl Stearate, Methylparaben, Dichlorobenzyl Alcohol.